


Tapping Massager with Far Infrared Heater

Model: 603F



INSTRUCTION MANUAL

INTRODUCTION

Thank you for purchasing the Tapping Massager with Far Infrared Heater. Combining the advantage of the Rubber Percussion Pad with Far Infrared Heat and the easy-to-use Rubber feel Handle, this is a perfect tool for soothing sore muscles and slimming body shapes, also improving body functions. The powerful percussion action works wonderful on your body. There are "High" and "Low" choices for percussion speeds; you could switch the Far Infrared Heater "On" or "Off" as you need. This percussion massager is going to give you a whole new experience of massaging.

With the high quality and the scientific design, the Rubber Percussion Head helps increasing the effectiveness in massaging. Also the Far Infrared Heater with 6-14uM Wavelength Far Infrared could bring you a great feeling of fast effectiveness of comfort result and far infrared ceramic heater works without consuming oxygen.

Please read this operating manual before using the Tapping Massager with Far Infrared Heater thoroughly for the correct methods of usage especially the Notes on Safety.

This manual should be kept available for future reference after you have read it thoroughly.

IMPORTANT SAFETY INSTRUCTIONS

WHEN USING AN ELECTRICAL APPLIANCE, BASIC PRECAUTIONS SHOULD ALWAYS BE FOLLOWED, INCLUDING THE FOLLOWING:

READ ALL INSTRUCTIONS BEFORE USING THIS MASSAGER AND SAVE THESE INSTRUCTIONS.

DANGER- To reduce the risk of electric shock

1. Always unplug the power cord from the electrical outlet immediately after using and before cleaning.

2. Never use or store the massager in a wet or moist condition; Never use it while taking a bath or shower.
3. If the massager has fallen into the water, unplug the power cord from the electrical outlet immediately; do not reach into the water to retrieve the massager right away.
4. Do not place or store the massager where it can fall or be pulled into a tub or sink. do not place in or drop into water or other liquid.

WARNING- To reduce the risk of burns, fire, electric shock, or injury to persons:

1. Never leave the massager unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
3. When children or the disables use this massager, close supervision by an adult is necessary.
4. The massager is just intended for household use. Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
5. Never operate the massager if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return your massager to a service center for examination and repair.
6. To reduce the risk of electric shock. The massager has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not change the plug in any way.
7. Do not carry the massager by supply cord as a handle.
8. Keep the cord away from heated surfaces.
9. Never operate the appliance with the air openings blocked.

Keep the air openings free of lint, hair, and the like. Never operate on a soft surface such as bed or couch where the air openings may be blocked.

10. Never drop or insert any object into any opening.
11. Do not use outdoors.
12. To disconnect, turn all controls to the off position, then remove plug from outlet.
13. Do not sit on the massager or press your body weights against it.
14. Do not use on the sensitive skin and body parts (e.g. eyes, heart).
15. Beware of the hot surface around the percussion pad and the heater enclosure.
16. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.


BENEFITS OF THE TAPPING MASSAGER WITH FAR INFRARED HEATER

- * The Far Infrared heating an area of the body produces specific results, the percussion massage goes deep into different layers of muscles, relieving tension, tones muscles and reduces stress which can increase the level of comfort and healing rate, leaving you rejuvenated and relaxed.
- * Using Far Infrared Heat therapy before exercising warms up your body for a better protection against injury before and after exercising. It stimulates blood circulation, preparing the muscles for action. It is useful in relieving pain, reducing swelling, and increasing mobility in soft tissues and muscles.
- * When heat is applied to a body part, the surface temperature rises, so does the temperature of the blood in the tissues just below the skin. The blood vessels enlarge to help speed up the removal of warmed blood from the area. This allows cooler blood to arrive, thus prevents from overheating, as a result of this vasodilatation. There is greater flow of blood through the heated area that brings in more oxygen and moves out more wastes.
- * A carefully balanced weight distribution system ensures that the

percussion motion is concentrated at the nodes. You receive the maximum massage effects with minimum efforts. A rocking action combines with the scientifically designed contoured nodes, delivering a penetrating percussion deep into your muscles. This results in a relaxing and invigorating feeling.

DESCRIPTION OF CONTROL AND COMPONENTS

- 1) Balance Wing for helping holding and pressing
- 2) Fans ventilation holes
- 3) "Infrared Heater" On/Off switch
- 4) "Speed" High/Off/Low switch
- 5) Rubber feel handle
- 6) Power cord
- 7) Holes for far infrared warm air
- 8) Rubber percussion head




OPERATION


- (1) Before using please make sure the switches are located at OFF and plug the power cord into any AC outlet.
- (2) Hold the rubber feel handle firmly in your hand and slide the "Speed" switch to "High" or "Low" to turn on the massager.
- (3) You may use the massager in standing, sitting, or lying down positions. Apply it on your neck, shoulders, back, waist thigh, calves or soles. (as figures)
- (4) Slide the "Heater" switch to "On" if there is need for Far Infrared Heat Therapy.




Shoulder Massage




Back Massage




Waist Massage



Thigh Massage



Calf Massage



Sole Massage

Note: Make sure that the fans ventilation holes are not blocked, avoid any small objects from falling into the fans ventilation holes. It is normal that the surface around the percussion pad and the fans area become warm after turning "On" the "Heater" switch.

- (5) While using the massager press it lightly on your body, do not press it too hard, as the device weight itself will help to give your body natural pressure of massaging.
- (6) During the usage you could move your body around and also switch the massager to different speed.
- (7) While finish using, please slide the switches "Off" and unplug the power cord.

MAINTENANCE

- * Use a dry cloth when cleaning the massager. Never use a wet cloth, use a cloth lightly dampened with water if desired.
- * Do not use alcohol, benzene, thinner or other harsh chemicals to clean the device.
- * Do not drop the massager.
- * Do not disassemble any parts and any other servicing should be performed by an authorized service representative or that has no user serviceable parts.
- * Avoid using the massager outdoor and do not leave it in direct sunlight or high humidity for a long period of time.
- * It is advisable that the power cord of the massager be coiled and the massager stored in a dust and moisture free environment, to prevent from affecting the functions.
- * Be sure to unplug the power cord after using and before cleaning.
- * Specific instructions for the people method of storage of the cord, the total appliance, and the like, when the appliance is not in use, and for care of the cord while in use (for example, for a hand-supported product, untwisting).